



THE MEN'S COLLECTION

SHAPING YOUR TRANSFORMATION



100% Manufactured in the USA
CELEBRATING 20+ YEARS OF SERVICE AND INNOVATION

WHY DO YOU NEED RAINEY RECOVERY WEAR?

Rainey Recovery Wear Compression Garments reduce the bruising and swelling during the early stages of recovery and promotes faster healing. Our garments deliver smoother results to the new contours of the patient's body.







1-2-3 POST-OP RECOVERY SYSTEM

As a patient progresses through each recovery step, they will need to transition to a different garment as the body's swelling is reduced. Each recovery step offers a garment with a different design function and feature. Rainey Recovery Wear's 1-2-3 Post-Op Recovery System ensures consistent compression throughout the entire healing process.

	STEP 1 SUPPORT	STEP 2 SCULPT	STEP 3 SLIM
WHEN TO WEAR	1 hour after surgery	1 week after surgery	4 weeks after surgery
DURATION OF WEAR	1-2 weeks	1-6 weeks	4-12 weeks
HOURS OF WEAR EACH DAY	24 hours	24 hours	12 hours
DESIGN FEATURE	DUAL-SIDE ZIPPER Simple donning and accessibility	ADJUSTABLE FRONT HOOK & EYE Adjust to tailored fit and delivers more compression	EASY PULL-ON Sized tighter to provide a long-term smoother slimmer result

BETTER FABRIC

Rainey Recovery Wear's Medical Grade Compression Fabric is engineered to deliver firmer compression than competitors. Our fabric is breathable, lightweight, and features a 4-way stretch, making it feel like a second skin!

-  Antimicrobial
-  Moisture Wicking
-  Control 4-Way Stretch
-  45% Nylon | 55% Lycra

Machine Washable on Cold and Air Dry



100% Manufactured in the USA



1-2-3 POST-OP RECOVERY SYSTEM

1

SUPPORT

A patient's body is swollen and sore and the incisions are seeping fluid. This garment is designed with Dual Side Zippers for easy donning and accessibility to maintain incision drainage pads and delivers the right amount of compression for immediate post-op recovery.

2

SCULPT

As a patient's body swelling continues to reduce, a patient needs to be able to adjust compression levels to their newly contoured body. The Front Hook & Eye design adjusts smaller to provide additional compression to create smoother results as a patient continues to heal.

3

SLIM

As a patient begins to engage in normal activity again, this garment provides additional compression to help maintain long-term results. This garment is easy to pull on and delivers more smoothing, slimming compression while maintaining a patient's smaller contoured body.





MEN'S FULL BODY

Rainey Recovery Wear Men's Full Body garments are excellent for post-op recovery from surgical procedures such as Male Breast Reduction, Gynecomastia, Body Contouring, Liposuction, Abdominal Liposuction, Abdominoplasty and other procedures focusing on the upper back, chest, upper and lower abdomen, flanks, or hips. Coverage varies.

STYLES

- Step 1: Support
- Step 2: Sculpt
- Step 3: Slim
- Short Sleeve

All styles are available in Boxer Shorts, Mid Thigh, Knee Length, Calf Length, and Ankle Length.

STEP 1: SUPPORT

Designed with Dual Side Zippers for easy donning and accessibility to maintain incision drainage pads and delivers the right amount of compression for immediate post-op recovery. Can be worn directly after surgery.



BOXER SHORTS

ITEM NO.
MF1-BX



MID THIGH

ITEM NO.
MF1-MT



KNEE LENGTH

ITEM NO.
MF1-KL



CALF LENGTH

ITEM NO.
MF1-CL


OPTIONS

Available Colors




Nude Black

Available Sizes: XS-6XL
Call for custom sizes.



ADJUSTABLE HOOK & EYE

Located on both shoulders ensures ideal fit and comfort.



DUAL-SIDE ZIPPERS

Easy donning and convenient access for both care of incisions and maintenance of post-op drainage pads.

NEW SUPPORT POUCH DESIGN

with Horizontal Fly
+ Full Open Crotch

Drain Site available upon request.

INSTRUCTIONS



To be worn 1 to 2 weeks



To be worn 24/7

Healing times may vary based on general health and medical maintenance.

ANKLE LENGTH

ITEM NO.
MF1-AL

STEP 2: SCULPT

As a patient's body swelling continues to reduce, a patient needs to be able to adjust compression levels to their newly contoured body. The Front Hook & Eye design adjusts smaller to provide additional compression to create smoother results as a patient continues to heal. Step 2 garments help achieve optimal results. Can be worn 1 week after surgery.



BOXER SHORTS

ITEM NO.
MF2-BX



MID THIGH

ITEM NO.
MF2-MT



KNEE LENGTH

ITEM NO.
MF2-KL



CALF LENGTH

ITEM NO.
MF2-CL


OPTIONS

Available Colors




Nude Black

Available Sizes: XS-6XL
Call for custom sizes.



ADJUSTABLE HOOK & EYE

Located on both shoulders ensures ideal fit and comfort.



FRONT HOOK & EYE CLOSURE

Offers 3 levels of compression that adjust progressively tighter to deliver firmer compression and support exactly where you need most!

NEW SUPPORT POUCH DESIGN

with Horizontal Fly
+ Full Open Crotch

ANKLE LENGTH

ITEM NO.
MF2-AL

INSTRUCTIONS



To be worn 1 to 6 weeks



To be worn 24/7

Healing times may vary based on general health and medical maintenance.

STEP 3: SLIM

As a patient begins to engage in normal activity again, this garment provides additional compression to help maintain long-term results. The easy pull on design delivers more smoothing, slimming compression while maintaining a patient's smaller contoured body.



BOXER SHORTS

ITEM NO.
MF3-BX



MID THIGH

ITEM NO.
MF3-MT



KNEE LENGTH

ITEM NO.
MF3-KL



CALF LENGTH

ITEM NO.
MF3-CL


OPTIONS

Available Colors



Nude Black

Available Sizes: XS-6XL
Call for custom sizes.



ADJUSTABLE HOOK & EYE

Located on both shoulders ensures ideal fit and comfort.


EASY PULL-ON DESIGN


Sized tighter to provide a long-term smoother, slimmer result.

NEW SUPPORT POUCH DESIGN

with Horizontal Fly
+ Full Open Crotch

INSTRUCTIONS

 To be worn 4 to 12 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.

ANKLE LENGTH

ITEM NO.
MF3-AL





ACTIVELY COLLABORATE WITH DOCTORS AROUND THE WORLD

Throughout Rainey Recovery Wear's history, our foundation has been set in research, science, and collaboration with practicing surgeons. Our founder and CEO, Russell Rainey has worked with surgeons around the globe to develop new, innovative garments for the procedures they are performing.

FULL BODY SHORT SLEEVE

The Men's Full Body Short Sleeve offers three levels of compression that adjust tighter as swelling is reduced to deliver firmer and evenly distributed compression, which helps skin re-attach smoothly to muscle tissue as a patient transforms. Can be worn directly after surgery.



BOXER SHORTS

ITEM NO.
MFS-BS

MID THIGH

ITEM NO.
MFS-MT

KNEE LENGTH

ITEM NO.
MFS-KL

CALF LENGTH

ITEM NO.
MFS-CL

OPTIONS



Available Colors



Nude Black

Available Sizes: XS-6XL

Call for custom sizes.



ADJUSTABLE FRONT HOOK & EYE CLOSURE

Offers 3 levels of compression that adjust progressively tighter to deliver firmer compression and support exactly where a patient needs it most.

NEW SUPPORT POUCH DESIGN


with Horizontal Fly
+ Full Open Crotch


Drain Site available upon request.

ANKLE LENGTH

ITEM NO.
MFS-AL

INSTRUCTIONS

 To be worn 1 to 6 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.





MEN'S UPPER BODY


Rainey Recovery Wear Men's Upper Body garments are excellent for post-op recovery from surgical procedures such as Abdominal Liposuction, Male Breast Reduction, Gynecomastia, and other procedures focusing on the chest, upper back, shoulders, lower abdomen, or flanks. Coverage varies.

STYLES

- Vest
- Chest Vest
- Vest Short Sleeve


MEN'S UPPER BODY

The Men's Upper Body garment offers three levels of compression that adjust tighter as swelling is reduced to deliver firmer and evenly distributed compression, which helps skin re-attach smoothly to muscle tissue as a patient transforms.



ADJUSTABLE HOOK & EYE

2-inch wide adjustable hook & eye fasteners at the shoulder ensure both ideal fit and comfort.




FRONT HOOK & EYE CLOSURE

Offers 3 levels of compression that adjusts progressively tighter to deliver firmer compression and support.

INSTRUCTIONS

 To be worn 1 to 4 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.



VEST SHORT SLEEVE

ITEM NO.
MVS


OPTIONS

Available Colors



Nude Black White (upon request)


Available Sizes: S-6XL
Call for custom sizes.




FRONT HOOK & EYE CLOSURE

Offers 3 levels of compression that adjust progressively tighter to deliver firmer compression and support.

INSTRUCTIONS

 To be worn 1 to 4 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.



UNISEX CONSULTATION ROBE

This elegant unisex Consultation Robe will make patients feel pampered and confident about their decision to move forward with their transformation. Rainey Recovery Wear's robe offers privacy and dignity for patients while providing easy access for the attending physician during consultation.



ITEM NO.
UR



OPTIONS

Available Color



Black

Available Sizes: S/M and L/XL
Call for custom sizes.

CARE

Rainey Recovery Wear Consultation Robe is easy to maintain with no dry-cleaning required. Just toss it into the washing machine (cold-water setting) and remove promptly from the dryer for wrinkle-free results.

FABRIC

100% polyester. Guaranteed not to fade or lose its shape.





FACIAL WRAPS

Rainey Recovery Wear Facial Wraps are excellent for post-op recovery from surgical procedures such as face lift, facial liposuction, chin augmentation, or cheek augmentation. Coverage varies.

STYLES

- Facial Wrap Non-Slip Dome Open Ears
- Facial Wrap Covered Ears
- Facial Wrap Non-Slip Dome Covered Ears
- Full Facial Wrap Covered Ears
- Facial Band Covered Ears
- Facial Wrap Open Ears

FACIAL WRAP NON SLIP DOME OPEN EARS

The Facial Wrap Non Slip Dome Open Ears delivers firm support and healing compression to reduce bruising and swelling.



DESIGN

Exclusive Non Slip Dome design keeps facial wrap in place during daily activities and while sleeping.

OPTIONS

Available Colors



Nude Black White

Available Sizes:


S, M, L, XL

INSTRUCTIONS

Wear as extended compression wear.

FACIAL WRAP COVERED EARS

The Facial Wrap Covered Ears is designed to adjust at both the head and neck to deliver the perfect level of firm support and therapeutic compression needed to reduce bruising and swelling, and ensure faster healing.



DESIGN

Adjustable Velcro Closures at head and back.

OPTIONS

Available Colors




Nude Black White

Available Sizes:

One size fits all

INSTRUCTIONS

 To be worn 1 to 3 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.

FACIAL WRAP NON SLIP DOME COVERED EARS

The Facial Wrap Non Slip Dome Covered Ears can adjust at the neck to deliver firm support and healing compression to reduce bruising and swelling.



DESIGN

Exclusive Non Slip Dome design keeps facial wrap in place. Adjustable Velcro Closures at the neck.

OPTIONS

Available Colors




Nude Black White

Available Sizes:

S, M, L, XL

INSTRUCTIONS

 To be worn 1 to 3 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.

FULL FACIAL WRAP COVERED EARS

The Full Facial Wrap Covered Ears can adjust at the head and neck to deliver the firm support and compression needed to reduce bruising and swelling for faster healing.



DESIGN

Adjustable Velcro Closures
at head and back.

OPTIONS

Available Colors




Nude Black White

Available Sizes:

S, M, L, XL

INSTRUCTIONS

 To be worn 1 to 3 weeks

 To be worn 24/7

Healing times may vary based
on general health and medical
maintenance.

FACIAL BAND COVERED EARS

The Facial Band Covered Ears can adjust at the crown of the head to deliver firm support and healing compression to help reduce bruising and swelling.



DESIGN

Adjustable Velcro Closures at head.

OPTIONS

Available Colors




Nude Black White

Available Sizes:

One size fits all

INSTRUCTIONS

 To be worn 1 to 3 weeks

 To be worn 24/7

Healing times may vary based on general health and medical maintenance.

FACIAL WRAP OPEN EARS

The Facial Wrap Open Ears can adjust at the head and neck to deliver the firm support and compression needed to reduce bruising and swelling for faster healing.



ITEM NO.
UF85



BACK VIEW

DESIGN

Adjustable Velcro Closures at head and neck.

OPTIONS

Available Colors



Nude

Black

White

Available Sizes:

One size fits all

INSTRUCTIONS



To be worn 1 to 3 weeks



To be worn 24/7

Healing times may vary based on general health and medical maintenance.

FACIAL WRAP NON SLIP DOME OPEN EARS

With Adjustable Velcro Closures

The Facial Wrap Non Slip Dome Open Ears can adjust at the neck to deliver firm support and healing compression to reduce bruising and swelling.



DESIGN

Adjustable Velcro Closures
at neck

OPTIONS

Available Colors



Nude Black White

Available Sizes:

S, M, L, XL

INSTRUCTIONS



To be worn 1 to 3 weeks



To be worn 24/7

Healing times may vary based
on general health and medical
maintenance.







AB BINDER

Rainey Recovery Wear Ab Binders are excellent for post-op recovery from abdominal liposuction, abdominoplasty, tummy tuck and other procedures focusing on the abdomen or flanks.

STYLES

- 9-inch Ab Binder
- 12-inch Ab Binder

9-INCH AB BINDER

The 9-inch Ab Binder features plush non-roll elastic with adjustable Velcro closures at front for a comfortable fit and firm targeted compression.



DESIGN

Delivers form-fitted support to the abdominal area.

OPTIONS

Available Colors



White

Available Sizes:

S/L and XL/2XL

INSTRUCTIONS



To be worn 4 weeks



To be worn 24/7

Healing times may vary based on general health and medical maintenance.

12-INCH AB BINDER

The 12-inch Ab Binder features plush non-roll Velcro compatible elastic with adjustable Velcro closures at front for a comfortable compression.




DESIGN

Delivers form-fitted support to the abdominal area.

OPTIONS

Available Colors



White

Available Sizes:
S/L and XL/2XL

INSTRUCTIONS



To be worn 4 weeks



To be worn 24/7

Healing times may vary based on general health and medical maintenance.




THE BEST FIT FOR THE BEST RESULTS

We guarantee each Rainey Recovery Wear garment is constructed to last the entire post-recovery period or you'll receive 100% reimbursement.

SIZING CHART

For the best fit, pull measurement tape medium snug when taking measurement.

UNISEX SIZING CHART				
Facial	S	M	L	XL
	23" – 24"	24.25" – 25.25"	25.5" – 26.5"	26.75" – 27.75"
AB Binder	S/M/L	XL - 2X		
	24" – 37"	38" – 47"		
Consultation Wear	S/M	L/XL		
Robe	24" – 37"	38" – 47"		



MEN'S SIZING CHART										
	S	M	L	XL	1X	2X	3X	4X	5X	6X
UPPER BODY										
Chest	36" – 38.5"	39" – 41.5"	42" – 44.5"	45" – 47.5"	48" – 50.5"	51" – 54.5"	55" – 58.5"	59" – 62.5"	63" – 66.5"	67" – 70.5"
LOWER BODY										
Waist (at navel)	30" – 32.5"	33" – 35.5"	36" – 38.5"	39" – 41.5"	42" – 44.5"	45" – 48.5"	49" – 52.5"	53" – 56.5"	57" – 60.5"	61" – 64.5"
Hips (at widest point)	36" – 38"	39" – 41"	42" – 44"	45" – 47"	48" – 50"	51" – 54"	55" – 58"	59" – 62"	63" – 66"	67" – 70"
Thigh (mid point)	19.75" – 21"	21 – 22.25"	22.25 – 23.5"	23.5 – 24.75"	24.75" – 26"	26" – 27.75"	27.75 – 29.5"	29.5 – 31.25"	31.25" – 33"	33" – 34.75"



ORDER ONLINE TODAY!

Same Day Delivery Available in the Atlanta Metro Area
www.raineywear.com

CONTACT INFO

Phone: 1-888-772-4639
Fax: 1-888-772-4630
Monday - Friday
8:30 AM - 5:30 PM EST



100% Manufactured in the USA